Cenefom

Cenefom Corp.

1F., No.50-5, Keyan Rd., Hsinchu Science Park, Zhunan Township, Miaoli County 35053, Taiwan

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1. **Product Description:**

The ENT sponge is provided sterile and for single use only. The device consists of PVA (polyvinyl alcohol) material which is capable of absorbing fluid several times of its own weight.

Product	Applicable surgery or treatment
Nasal Packing Standard	Septoplasty, turbinectomy, rhinoplasty surgery
Nasal Packing Epistaxis Type Nasal Packing Nosebleed Type	Nosebleed treatment
Nasal Packing Sinus Type	Sinus surgery
Nasal Packing Anatomical Type	Nasal surgery to fill the nasal cavity
Ear Pack	Otic surgery
Ear Wick Ear Wick / Fenestrated	Otitis externa related surgery

2. **Intended Use:**

The ENT sponge is designed for placing inside the nasal cavities or ear canals as a tamponade or hemostasis material to control oozing or persistent bleeding.

3. **Intended Therapeutic Use and Indications:**

The ENT sponge is mainly used for hemostasis for ear canal (meatus), otitis chambers, bleeding within nasal cavity. Bleeding can be stopped effectively after insertion of ENT sponge to provide a tamponade effect. It is intended to be used either after performing septoplasty, sinus surgeries, other nasal cavity surgeries and otology operations of the external auditory canal operations such as canalplasty or myringoplasty. An appropriate length and specification of polymer packing (or also containing a silicone vent tube) can be placed for nasal and otology surgeries at the bottom of the meatus or nasal cavity (single side or bilaterally).

4. **Intended Patient Population:**

Male and female older than 7 years.

5. **Intended User Profile:**

- 5.1. Only licensed clinicians are allowed to perform placement or removal of the ENT sponge to prevent patient harm due to incorrect operation. If found the product loss of function, please stop using it immediately because it is impossible to be sure that the foam can effectively achieve its intended use.
- 5.2. The clinicians shall select appropriate size of the ENT sponge to place into nasal cavity or ear canal. It is recommended to use one ENT sponge on each side of the nose or ear. Place more than one ENT sponges on each side of the nose or ear should be considered by clinicians.
- 5.3. The ENT sponge shall be removal by clinicians within 72 hours after insertion.
- 5.4. The ENT sponge to sale by, or on the order of, a physician.

6. **Contraindication:**

- 6.1. Absolute Contraindication:
Patient allergy to PVA (polyvinyl alcohol), nylon, glycerin or silicon materials.
- 6.2. Relative Contraindication:
Pregnant woman or the patients who have had any adverse reactions by using PVA, nylon, glycerin or silicon related materials.

7. **Side Effects:**

- 7.1. An improper introduction and removal of the tamponade can lead to pain, bleeding and/or traumatization of the mucous membranes of the cavity inserted.
- 7.2. If the ENT sponge is placed more than 72 hours, it might lead to local infection or systemic infection such as toxic shock syndrome.
- 7.3. No other side effects are known if applied properly.

8. **Precaution:**

- 8.1. Do NOT re-sterilise or re-use. Devices are supplied sterile by irradiation ready for use.
- 8.2. Do not use if the packaging has been opened or damaged. Otherwise it may cause inflammation or systemic infection.
- 8.3. Do not reuse, as it may result in compromised device performance and risk of cross contamination. It might lead to local infection or systemic infection such as toxic shock syndrome.
- 8.4. The ENT sponge is not suitable to be moistened with an alcoholic solution as the structure of the sponge will be changed and resulting in loss of function, when the Product loss of function, it is impossible to be sure that the foam can effectively achieve its intended use.
- 8.5. The ENT sponge is not suggested to use on people who may be allergic to PVA sponge.
- 8.6. Please note that trimming the foam may cause risks, such as displacement and ingestion of a nasal packing may occur.
- 8.7. Due to the unique production process of the tubing products, there is a chance that the following situations may occur:
As the product is stored in vacuum packaging, once the product is opened and moistened and has expanded, there is the possibility that the two ends of the silicon tube may not expand automatically. To avoid this condition, please carefully press and adjust the ends of the silicon tube before moistening the foam.
Please note that the above condition will not affect the original functionality of the product.
- 8.8. The product comes into contact with bodily fluids, which can be contaminated. Dispose of used sponges according to local medical regulations to prevent the risk of environmental contamination or reuse.
- 8.9. Securing the string to the patient's cheek will prevent the sponge moving backwards, thus eliminating the possibility of the Nasal Packing being imbibed further and causing an obstruction.
- 8.10. Nasal Packing with tube type with a silicone tube that allows the patient to breathe easily. If finding that the silicone tube may become loose during use, please stop using it immediately because the silicone tube may get stuck in the trachea to cause harm.
- 8.11. Because children under the age of 7 are not fully developed, it is impossible to be sure that the foam can effectively achieve its intended use.
- 8.12. Normally, the edges of PVA foam are smooth. Do not use if it is not smooth to avoid causing pain to the surgical site and unintentional wounds.

9. **Use Procedure:**

Please note that the PVA sponge will only expand properly by aqueous solutions moisten (do not use oil or fat-containing solutions/drugs to moisten).

Step 1. Prior to placing the ENT sponge

- a. Select appropriate size of the ENT sponge, otherwise cannot be stopped effectively after insertion of ENT sponge to provide a tamponade effect.
- b. Placement into nasal cavity or ear canal.
(Caution:
 - A. The instruments of contacting sponge shall be clean and disinfect according the medical-related requirements.
 - B. When using sponge with string, the string must be on the outside.
 - C. The tied string position toward the D-shaped protrusion refer to the dotted circle of below figure. The placement direction must be correct, otherwise cannot be stopped effectively after insertion of ENT sponge to provide a tamponade effect.)

Step 2. Removal of the ENT sponge

- a. Prior to attempting retrieval of the ENT sponge, the clinician may apply normal saline and wait for few, no more than 10 minutes until the ENT sponge is thoroughly moistened.
- b. Using a tweezers to gently pull the ENT sponge. If significant resistance is encountered, repeat step 2.a and pull again. Care must be taken not to damage surrounding mucosa if excessive adhesion occurs.
(Caution: The string serves for fixation only, it must not be used for removing the sponge, otherwise it may cause harm.)
- c. Repeat step 2.b until the ENT sponge is completely retrieved from the place inserted.

10. Warnings:

- 10.1. Apply any medication onto the ENT sponge for the purpose of hemostasis or anti-adhesion shall be directed or performed by a physician.
- 10.2. After applying the product more than three days, the evaluation by physician is suggested for further usage.
- 10.3. If the patient experiences irritation, infection after insertion of a ENT sponge, the ENT sponge should be removed by the physician immediately.
- 10.4. Rarely, displacement and ingestion of a nasal packing may occur. PVA surgical sponge is non-toxic and is considered safe to pass once it is completely in the alimentary tract. The patient should be monitored in the normal manner for the ingestion of any non-toxic material. Clinical confirmation that packing material has not been aspirated or is not in danger of being aspirated should be assessed. If aspiration should occur, prompt medical intervention is imperative with close patient monitoring and airway supporting pending the removal of the material.
- 10.5. In rare instances the physiochemical conditions associated with nasal or meatus or otitis chambers surgery, both with and without ENT sponge, may present a risk of toxic shock syndrome (TSS). Warning signs of TSS include sudden fever, vomiting, diarrhea, dizziness, fainting (or near fainting when standing up), and/or a rash that looks like a sunburn the treatment with ENT sponge has to be interrupted immediately. ENT sponge exhibit no antimicrobial properties, they are not bacteriostatic towards pre-existing infections, nor do they prevent the occurrence of new infections. In the case of pre-existing infections it is recommended that appropriate treatment be instituted.
- 10.6. When removing nasal packing, make sure all parts are completely removed to avoid risk of aspiration to the patient. If aspiration should occur, prompt medical intervention is imperative with close patient monitoring and airway supporting pending the removal of the material.

11. Preservation & Handling:

- 11.1. The product should be preserved in a dry, well-ventilated place and in sealed conditions.
- 11.2. Sterility and performance are not guaranteed on devices that have passed their Use By date or if the packaging or seals have been damaged. If label including use by date is illegible, do not use.
- 11.3. Please storage between 0°C ~ 37°C to prevent material deterioration.

12. Incident Reporting:

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

13. Safety:

The ENT sponge has been tested and shown to be sterile, non-pyrogenic, non-cytotoxic, non-irritating and non-sensitization.

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