# Cenefom

Cenefom Corp.

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# EN English

## 1. **Product Description:**

The Ophthalmic Sponge is a sterile, single use device. The OP-ES series consisted of a spear is made of PVA (polyvinyl alcohol) material which can absorb aqueous fluids several times of its mass instantaneously upon contact with the fluids and a handle. The OP-EW, OP-CS, OP-LR, OP-LRD and OP-IW series are made of PVA (polyvinyl alcohol) material which can absorb aqueous fluids several times of its mass instantaneously upon contact with the fluids. They have a cumulative transient use of less than 60 minutes on the conjunctiva (mucous membrane of the eye).

Product Type	Applicable surgery or treatment
PVA Eye Spear or PVA Sponge Point	Absorbing excess fluids from the operative site.
Eye Wick	Drainage excess fluids and prevent them from entering the surgical site.
LASIK Ring	Absorbing and prevents the excess fluids and blood from encroaching into the surgical site.
LASIK Eye Drain	Absorbing excess fluids and protection of flap and hinge.
Corneal Light Shield	Keeping the cornea moist, protect the retina from potentially harmful light during ophthalmic procedures.
Wipe	Mopping fluids.

## 2. **Intended Use:**

The Ophthalmic Sponge is indicated for using in ophthalmic microsurgical procedures for absorbing or mopping fluids.

## 3. **Intended Therapeutic Use and Indications:**

The Ophthalmic Sponge is mainly used for absorbing aqueous fluids in ophthalmic or microscopic surgical procedures.

## 4. **Intended Patient Population:**

Male and female older than 7 years.

## 5. **Intended User Profile:**

- 5.1. The ophthalmic sponge should be used on the order of a physician.
- 5.2. The ophthalmic sponge to sale by, or on the order of, a physician.
- 5.3. Only licensed clinicians are allowed to perform use of the Ophthalmic Sponge to prevent patient harm due to incorrect operation. If found the product loss of function, please stop using it immediately because it is impossible to be sure that the foam can effectively achieve its intended use.

## 6. **Contraindications**

- 6.1. Absolute Contraindication:  
Patient allergy to PVA (polyvinyl alcohol), PP(Polypropylene) or glycerin materials.
- 6.2. Relative Contraindication:  
Pregnant woman or the patients who have had any adverse reactions by using PVA (polyvinyl alcohol), PP(Polypropylene) or glycerin related materials.

## 7. **Side Effects:**

- 7.1. Irritation if any loose particles are not rinsed away from the treatment site.
- 7.2. Foreign body granuloma if any loose particles are not rinsed away from the treatment site.
- 7.3. Bleb inflammation if any loose particles contaminate in or around a bleb.

## 8. **Precaution:**

- 8.1. Do NOT re-sterilise or re-use. Devices are supplied sterile by irradiation ready for use.
- 8.2. Do not use if the packaging has been opened or damaged, otherwise, it may cause inflammation.

- 8.3. Do not reuse, as it may result in compromised device performance and risk of cross contamination. It might lead to local infection or systemic infection.
- 8.4. The ophthalmic sponge is not suitable to be moistened with an alcoholic solution as the structure of the sponge will be changed and resulting in loss of function, when the Product loss of function, it is impossible to be sure that the foam can effectively achieve its intended use.
- 8.5. The ophthalmic sponge is not suggested to use on people who may be allergic to PVA sponge.
- 8.6. The product comes into contact with bodily fluids, which can be contaminated. Dispose of used sponges according to local medical regulations to prevent the risk of environmental contamination or reuse.
- 8.7. Because children under the age of 7 are not fully developed, it is impossible to be sure that the foam can effectively achieve its intended use. Please ask a doctor for careful evaluation.
- 8.8. PVA sponge may not absorb water as fast as expected due to some reasons. You can increase the temperature of the wetting liquid to 35~37 °C, or continue to soak the sponge in the wetting liquid until the foam is completely expanded.
- 8.9. Please choose the appropriate product according to the patient's condition, and make sure the foam is fully soaked to avoid causing harm to the patient.
- 8.10. If the product is damaged when the package is opened, do not use it, otherwise it may cause harm.

## **9. Use Procedure:**

### **9.1. PVA Eye Spear:**

#### **9.1.1. Opening The Product Package**

- 9.1.1.1. The product package contains two foiled pouches sealed in an envelope. The envelope has a Tyvek top layer and a clear plastic bottom layer.
- 9.1.1.2. To open the package, separate the clear plastic layer from the Tyvek layer and gently pull open the envelope.
- 9.1.1.3. Empty the two foiled pouches onto a sterile field.  
Caution: Maintain sterile technique when contacting the inner foiled pouches.

#### **9.1.2. Opening The Sterile Foiled Pouches**

- 9.1.2.1. Each foiled pouch contains five (5) PVA SPEAR.
- 9.1.2.2. While maintaining sterile technique, tear open the foiled pouch.
- 9.1.2.3. Remove spears and place in the sterile field.  
Caution: In order to avoid inadvertent contact of the spear with aqueous fluids, keep the spears in a dry place within the sterile field.

#### **9.1.3. Using The PVA Spear**

- 9.1.3.1. By holding the plastic handle, gently touch the tip of the spear at the site where excess fluid is to be removed.
- 9.1.3.2. Hold still for approximately 3 seconds to allow for complete swollen of the PVA material.  
Caution: Do not squeeze or press the PVA spear against the surface of the surgical sites. Might cause harm to the surgical site.
- 9.1.3.3. Dispose of the used spear into a biological waste container.
- 9.1.3.4. When necessary, use additional spears to absorb fluid until all excess fluid is satisfactorily removed.

### **9.2. Eye Wick:**

- 9.2.1. Moisten wick(s) with sterile fluid and insert around eye site (without pinching).
- 9.2.2. Irrigate normally.

### **9.3. LASIK Eye Drain & LASIK Ring:**

- 9.3.1. Remove LASIK Eye Drain/LASIK Ring and place in the sterile field. Moisten LASIK Eye Drain/LASIK Ring if needed.
- 9.3.2. Use the LASIK Eye Drain/LASIK Ring at the surgical process site where excess fluid is to be removed.

### **9.4. Corneal Light Shield:**

- 9.4.1. Remove Corneal Light Shield and place in the sterile field. Moisten Corneal Light Shield if needed.
- 9.4.2. Place the Corneal Light Shield onto the cornea to prevent light illumination and maintain cornea moist throughout surgery.

### **9.5. Wipe:**

- 9.5.1. Moisten wipe(s) with sterile fluid and mop up fluids during surgical procedures.  
Caution: Do not clean or removal of residue or other materials in and on the surgical site.

## **10. Warnings:**

- 10.1. Do not attempt to use the sponge for cleaning of the surgical sites or removing debris, such as vitreous

- material that might be generated during the surgery. Might cause harm to the surgical site.
- 10.2. The ophthalmic sponge is not for use as an applicator for applying medication to the surgical sites. Also, the ophthalmic sponge shall not be used as a hemostasis tool for compressing the spear against the bleeding site. Might cause harm to the surgical site.
  - 10.3. If any signs of irritation or infection to the eyeball or cornea are observed, do not use the ophthalmic sponge and consult a physician immediately.

**11. Preservation & Handling:**

- 11.1. The product should be preserved in a dry, well-ventilated place and in sealed conditions.
- 11.2. Sterility and performance are not guaranteed on devices that have passed their Use By date or if the packaging or seals have been damaged. If label including use by date is illegible, do not use.
- 11.3. Please storage between 0°C~ 37°C to prevent material deterioration.

**12. Incident Reporting:**

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

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